


1.5 METER WADSWORTH SPECTROGRAPH




A COMPACT, LOW COST, SPECTROGRAPH
OFFERING EXCELLENT PERFORMANCE,
VERSATILITY, AND OPERATIONAL ECONOMY



Since 1947 the Jarrell-Ash Company has offered a high quality, low cost 1.5 Meter Wadsworth Spectrograph. Its sound design and workmanship are time-tested in laboratories which represent such divergent fields as industry, medicine, agriculture, research, and the academic.

The Jarrell-Ash Company also manufactures a 3.4 Meter Ebert, Convertible Spectrograph. This large grating instrument provides more versatility with the best possible resolution of any commercial instrument.


The 1.5 Meter Wadsworth Spectrograph meets the need for a simple, compact, low cost, fixed position instrument of professional quality which is not a compromise, but will meet most laboratory requirements. One of the three models offered herein will fill these requirements.



FEATURES

FULL 20 INCHES OF QUALITY STIGMATIC SPECTRA

Unlike other concave grating mounts, the Wadsworth produces stigmatic spectra by using a concave mirror to illuminate the grating. The spectrum, therefore, has a point-to-point correspondence to the source, producing sharply defined lines having a high ratio of signal-to-background.




A neutral density filter or stepped sector can be used at the slit for wavelength calibration of the entire range of the instrument.

A "fishtail" can be used at the slit to set the height of the spectrum.

A Hartmann diaphragm can be used at the slit to compare a series of exposures.


CAMERA RACKING



The complete 20" camera is racked by a precision gear system controlled from the slit end of the spectrograph.


As many as twenty sets of spectra .9 mm high can be obtained on one film. Each set is the image of the same slit position. Different positions of the slit are not used to obtain successive exposures.

100 FOOT CASSETTE



For convenience, the full 100' roll of film is placed in a cassette which feeds into the film holder. The entire assembly can be brought into the dark room and the 20" strip cut off. An extra assembly can be used if more than one film type is employed routinely.

STAINLESS STEEL PRECISION SLITS



Because the spectrum recorded is the image of the slit, it is imperative for quality spectra that a sharp evenly spaced slit be employed. Bilateral adjustable or fixed bayonet type precision ground stainless steel (not evaporated) slits are employed.

SPEED

The f/24 spectrograph is illuminated with a highly effective cross cylindrical system to fill the collimator completely. High efficiency precision gratings ruled by Jarrell-Ash for the most effectiveness at the wavelengths required yield working spectra in fractions of a second.

STURDY CONSTRUCTION


Steel U-beams of $\frac{1}{8}$ " stock make up the frame of the spectrograph. A 16-gauge aluminum case provides convenient access ports to the grating and mirror.

PORTABILITY

Total weight of the spectrograph is only 325 lbs. A dolly is available for convenience in moving to various instruments.

TYPICAL SPECTRA

Taken with Standard Model 78-000 using a 590 grooves/mm grating.



S. Q. Powder spectrum using D. C. Arc and 1-N Emulsion with set of two filters. Upper spectrum shows 5600-8000A — Lower spectrum shows 2500-4000A

Taken with Wide Angle Model 78-090 using a 1180 grooves/mm grating.




Upper spectrum — 316 stainless steel and lower spectrum — low alloy steel (cont.)

Taken with Wide Angle Model 78-090 using a 590 grooves/mm grating.




Portion of Copper spectrum in Neon filled hollow cathode lamp using D.C. source and SA-1 emulsion.

Portion of Copper spectrum in Neon filled hollow cathode lamp using D.C. source and SA-1 emulsion.



THE 1.5 METER WADSWORTH SPECTROGRAPH FOR RESEARCH

The requirement for a portable spectrograph of moderate speed with a broad wavelength coverage is often found in a research laboratory. The 1.5 Meter Wadsworth meets this requirement. Mounted on an adjustable dolly, the Wadsworth can be wheeled up to an experiment and quickly obtain the required spectra.



The Wadsworth used for plasma studies at Sperry Rand Corporation, Sudbury, Massachusetts.


Spectra of short-life species occurring in a shock tube can be recorded. The broad wavelength coverage guarantees that nothing will be missed. Since only one shot is often all that will be taken, the lower cost of film is an obvious advantage.

Other typical uses include combustion studies of flames, temperature measurements of plasma arcs, missile exhausts, incandescent materials. The Wadsworth is the ideal instrument for the varied requirements of a research operation.



Emission spectral studies of combustion in a rocket engine with a 1.5 Meter Wadsworth at NASA Lewis Research Center, Cleveland, Ohio.

Auxiliary apparatus such as the model 18-022 optical scanning condensing system enables portions of an extended source to be studied.





THE WADSWORTH 1.5 METER SPECTROGRAPH FOR CHEMICAL ANALYSIS

Routine multi-element analysis with the Wadsworth at the Eastern Stainless Steel Corporation, Baltimore, Maryland.




The Jarrell-Ash Wadsworth Spectrograph is currently used in many diverse fields as an analytical tool. It was chosen by these laboratories because it met their particular needs at the lowest cost. The versatility of the Wadsworth Spectrograph is demonstrated by the diversity of the industries represented, such as chemical (both organic and inorganic), biological, toxicological, geological, metal and semi conductor, academic, forensic, testing labs and even art museums.

The mode of exciting the sample is as varied as the laboratories using the spectrograph. It is most important to choose the correct source to go with your spectrograph. The Jarrell-Ash Varisource is a time tested source which provides a large number of parameters, convenient and safe operation. Types of excitation available are: DC Arc, High Voltage AC Spark, Unidirectional DC Arc (Half wave and full wave), Low Voltage DC Arc, Self-Sustaining Low Voltage AC Arc, and High Voltage AC Arc. Safety lock controls also enable the operator to choose the resistance capacitance and inductance required. For full particulars see the Varisource catalog.




LASER MICROPROBE SOURCE

The laser microprobe is a unique mode of sample excitation. The analysis of small spots or inclusions (50 microns) in conducting or non-conducting materials can be quickly performed with virtually no previous sample preparation. Metals, refractories, biological materials all can be analyzed. See the Laser Microprobe catalog for full particulars.



The Laser Microprobe/Wadsworth System used for ultra-micro analysis of metallics in biological tissue at Stanford University, School of Medicine, Palo Alto, California.

ARC-SPARK STAND



A convenient, easy-to-use arc and spark stand enables pins, flats, and solution samples to be excited. The stand is designed so that only the electrode clamps project into the chamber. Removable G.S.G. silicone inserts are used to line the chamber. These features reduce the possibility of contamination and permit easy cleaning. All control knobs are located in logical, easy-to-use positions. The electrodes can be viewed during exposure through a polarizer window at 6x magnification. The electrode jaws open bilaterally. They are made of nickel-plated brass coated with Rhodium to minimize

any possibility of corrosion and contamination. The electrode clamps are internally cooled to withdraw any heat generated. Attachments for handling various forms of samples are available.

2. ROTATING DISC OR ROTATING PLATFORM ATTACHMENT

For solution samples, the Model 19-330 solution attachment is installed in the arc stand. Solutions are placed in a ceramic boat. A disc type of electrode dips into the solution. Rotation of the disc brings a fresh solution to the counter electrode at a rate dependent on the drive motor employed.


The rotating platform attachment Model 19-340A enables samples such as powders which are formed in disc shape to be rotated such that a fresh area is exposed to the counter electrode. This technique is also used when thin samples are encountered.

3. PETRY STAND

Flat samples of varying sizes and shapes are often encountered. The Model 19-310 Petrey stand is the conventional method for holding such samples. The counter electrode is placed in the lower clamp and a fixed gap spacer is used to assure repeatable spacing.

SPECIFICATIONS

	STANDARD MODEL	WIDE ANGLE MODELS	
CATALOG NUMBER	78-000	78-090	78-090
RULING OF GRATING			
GROOVES/mm	590	590	1180
WAVELENGTH RANGE A			
First Order	2100-7800	4200-9600	2100-4800
Second Order	2100-3900	2100-4800	2100-2400
Third Order	2100-2600	2100-3200	
RECIPROCAL LINEAR DISPERSION			
A/mm			
First Order	10.9	10.8	5.4
Second Order	5.45	5.4	2.7
Third Order	3.6	3.6	
THEORETICAL RESOLVING POWER			
First Order	33,700	33,700	67,000
GUARANTEED RESOLUTION	.1A	.1A	.1A
NET WEIGHT	295	315	315
SHIPPING WEIGHT	365	425	425



PARTIAL LIST OF USERS

AIR REDUCTION COMPANY, INC.
Murray Hill, New Jersey
AMERICAN CYANAMID COMPANY
Stamford, Connecticut
ATLANTIC REFINING COMPANY
Philadelphia, Pennsylvania
BROOKHAVEN MEMORIAL HOSPITAL
East Patchogue, New York
BROWN UNIVERSITY, DEPT. OF PHYSICS
Providence, R.I.
EDGERTON, GERMESHAUSEN AND GRIER
Boston, Massachusetts
CARBORUNDUM COMPANY
Kearse, New Jersey
MASER OPTICS CORPORATION
Cambridge, Massachusetts
CARRIER CORPORATION
Syracuse, New York
M.I.T. MINERAL ENGINEERING LAB.
National Lead Company
Hightstown, New Jersey
MUSEUM OF FINE ARTS
Boston, Massachusetts
EASTERN STAINLESS STEEL CORP.
Baltimore 24, Maryland
U.S. ARMY SIGNAL CORPS.
Fort Monmouth, New Jersey
NATIONAL RESEARCH CORPORATION
Cambridge, Massachusetts
GENERAL ELECTRIC COMPANY
Schenectady, New York
PRATT AND WHITNEY AIRCRAFT
Middletown, Connecticut
GENERAL ELECTRIC COMPANY
Syracuse, New York
RAYTHEON COMPANY
Waltham, Massachusetts
WILLIAM GOLD REFINING CO., INC.
Buffalo 14, New York
WESTINGHOUSE ELECTRIC CORP.
Springfield, Massachusetts
HARBISON-WALKER REFRACTORIES CO.
Pittsburgh, Pennsylvania
I.B.M. CORP.
Vestal, New York
PHELPS DODGE COPPER PRODUCTS
South Brunswick, New Jersey
JOHNS HOPKINS UNIVERSITY
Applied Physics Laboratory
OLMSTED AIR FORCE BASE
Middletown, Pennsylvania
LEHIGH TESTING LABORATORY
Wilmington, Delaware
A. S. McCREATH & SONS TESTING LAB.
Harrisburg, Pennsylvania

NAT. AERONAUTICS & NASA, LANGLEY
Hampton, Virginia
VANDERBILT UNIVERSITY
Department of Chemistry
STATE OF ALABAMA
Dept. of Toxicology
and Criminal Investigation
Mobile, Alabama
UNIVERSITY OF VIRGINIA
Department of Physics
SHELL CHEMICAL COMPANY
Norco, Louisiana
PIPELINE MAINTENANCE COMPANY
Tulsa, Oklahoma
NORTH CAROLINA STATE
Board of Health
Raleigh, North Carolina
LOUISIANA STATE UNIVERSITY
New Orleans, Louisiana
MISSISSIPPI STATE UNIVERSITY
State College, Mississippi
TERNSTEDT DIVISION GMC
Detroit 9, Michigan
ARGONNE NATIONAL LABORATORY
TWIN CITY TESTING & ENGINEERING LAB.
St. Paul, Minnesota
ARGONNE CANCER RESEARCH HOSPITAL
Chicago, Illinois
STANDARD OIL OF OHIO
Cleveland, Ohio
PURDUE UNIVERSITY
ALCOA
St. Louis, Missouri
ENGLES SHEPPARD DIV. OF WYMAN-GORDAN
Harvey, Illinois
LITTON INDUSTRIES
Minneapolis, Minnesota
PENNSALT CHEMICAL CORPORATION
Wyandotte, Michigan
NATIONAL AERONAUTICS NASA, LEWIS
Cleveland, Ohio
MONARCH ALUMINUM MFG. CO.
Cleveland, Ohio
MISSISSIPPI LIME COMPANY
St. Genevieve, Missouri
MIDWEST RESEARCH INSTITUTE
Kansas City, Missouri
UNIVERSITY OF MICHIGAN
Pontiac, Michigan
WELLMAN BRONZE & ALUMINUM CO.
Cleveland, Ohio
SANDIA CORPORATION
Albuquerque, New Mexico

PHYSICIANS MEDICAL LABORATORY
Portland, Oregon
NORTH AMERICAN AVIATION
Chatsworth, California
AEROSPACE CORPORATION
El Segundo, California
ARIZONA STATE UNIVERSITY
THE MARTIN COMPANY
Englewood, Colorado
LOCKHEED AIRCRAFT CORPORATION
Palo Alto, California
EL PASO NATURAL GAS COMPANY
Farmington, New Mexico
STAUFFER CHEMICAL COMPANY
Henderson, Nevada
U.S. GEOLOGICAL SURVEY
Denver 25, Colorado
WESTINGHOUSE ELECTRIC CORP.
Pittsburgh 35, Pennsylvania
CALIFORNIA INSTITUTE OF TECHNOLOGY
NATIONAL BUREAU OF STANDARDS
Boulder, Colorado
DEPT. OF MINES & TECHNICAL SURVEYS
Ottawa, Ontario, Canada
HARVARD UNIVERSITY
HERCULES POWDER COMPANY
Bacchus, Utah
OPTICS TECHNOLOGY, INC.
Belmont, California
ESCUELA TECNICA SUPERIOR DE
INGENIEROS INDUSTRIALES
Bilbao, Spain
INSTITUT FUER GERICHTLICHE MEDIZIN
DER UNIVERSITAET ERLANGEN
Erlangen, Germany
MITSUBISHI METAL MINING CO., LTD.
Osaka, Japan
INSTITUTE OF MINERALOGY
Basel, Switzerland
MINERAL SEARCH OF AFRICA, LTD.
Northern Rhodesia, Africa
UNIVERSIDAD NACIONAL DEL LITORAL
Rosario, Argentina
UNIVERSITAET MUENCHEN, INSTITUT
FUER ALLGEMEINE UND
ANGEWANDTE PHYSIK
Muenchen 2, Germany
MONTECATINI S. p. A.
Scarlino Follonica, Italy
GREEK ATOMIC ENERGY COMMISSION
Athens, Greece
UNIVERSITY OF QUEENSLAND
Australia

NATIONAL AERONAUTICS & NASA, LEWIS



JARRELL-ASH COMPANY

590 LINCOLN STREET, WALTHAM, MASSACHUSETTS, 02154
TELEPHONE: 617-899-4300 CABLE: JACO-WALTHAM

SAN MATEO, CAL.
1789 Lexington Avenue
415 Fireside 1-2814

COSTA MESA, CAL.
2533 Fordham Drive
714 545-7350

ATLANTA, GA.
1430 Monroe Drive, N.E.
404 TRinity 6-9020

HOUSTON, TEXAS
Southwestern Engineering
& Equipment Co.
4101 San Jacinto Street
713 Jackson 8-5538

DALLAS, TEXAS
Southwestern Engineering
& Equipment Co.
Suite 105
6300 N. Central Expressway
214 EMerson 3-6501

FARMINGTON, MICH.
34361 Quaker Valley Road
313 476-5936

ELMHURST, ILL.
110 East Schiller Street
312 834-2878

WASHINGTON, D. C.
933 G Street, N. W.
202 628-3908

NEW BRUNSWICK, N. J.
347 Livingston Avenue
201 Kilmer 5-5305

PITTSBURGH 9, PA.
916 Greenhill Road
412 821-2990

JENKINTOWN, PA.
628 Pinetree Road
215 885-1055

CANADA (Except British Columbia) TECHNICAL
SERVICE LABORATORIES, 355 King Street, West,
Toronto 2B, Ontario, 416 EMpire 2-4248
(British Columbia served by San Mateo Office)

NIPPON JARRELL-ASH CO., LTD.
Kiyomachi-Sanjo-Sagaru
Nakagyo-ku
KYOTO, JAPAN

JARRELL-ASH (EUROPE) S. A.
Rue de la Jaluse, 6
LE LOCLE, SWITZERLAND